Palodent® Plus

Système de Matrices Sectorielles

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With the tines properly positioned interproximally, release the locking slide, and slowly release the forceps tension. When ring is securely seated, release the forceps completely and withdraw forceps.

Small parts under tension

- Be sure ring is securely seated in the forceps before introducing into the patient field.
- Be sure ring is securely seated apically before withdrawing forceps.

4. Cavity Preparation

Complete cavity preparation as usual. Once the preparation is complete, remove the ring or the metal shield of the WedgeGuard (if applicable). Use the pin tweezers to grasp the hole on the WedgeGuard shield, and slide shield laterally to remove. The wedge may remain in place if apical to the gingival cavsosurface margin. Properly dispose of the shield with contaminated sharps.

5. Matrix Selection

The 5.5mm Palodent® Plus contoured sectional matrix is recommended for most applications. The 4.5mm and 3.5mm are designed for use with adolescents and patients with poorly erupted posterior teeth. They are also ideal for the direct-indirect technique. The 6.5mm matrix is designed for larger cavities occlusogingivally. It augments a contoured matrix with a gingival apron and an adjustable continuation of the contour at the marginal ridge. The 7.5mm matrix is useful for patients with periodontal issues that require more coverage toward the gingival margin.

Technique Tip: Rings, Matrices, Wedges and WedgeGuards may be used with most popular matrix systems, including other sectionals. The Rings, Wedges and WedgeGuards may be used with most circumferential matrices. The benefits of the Ring, Wedge and WedgeGuard contours assist proper proximal contour.

6. Application of Matrix and Wedge

Grasp the matrix with the pin tweezers using the hole located on the occlusal tab. If the ring was placed properly before tooth preparation, then the matrix should slide easily into place adjacent to the space to be restored from the occlusal approach. Check to be sure gingival tissue is not trapped between band and tooth surface. If so, slightly lift and reposition band. Ensure matrix is apical to preparation margin with no gap. Technique Tip: If WedgeGuard wedge prevents matrix from fully seating gingivally, slide wedge out slightly, place matrix, then reposition wedge.

- Close the gingival margin by placing a contoured wedge or anatomic wooden wedge. Palodent® Plus wedges contain a placement tab with hole, which is securely grasped by the pin tweezers. This assures secure placement and minimizes the likelihood of dropping small accessories. Technique Tip: Palodent® Plus wedges' hollow design allows for easy placement of a second wedge from the opposing side (if needed) without dislodging the first. Technique Tip: Alternatively, the fork design allows placement of the Palodent® Plus retaining ring first if desired. The wedge can then be easily placed through the fork tine without dislodging the ring.

7. Replace the Palodent® Plus retaining ring in the proximal plane, engaging the matrix and the wedge, thus stabilizing it against tooth structure. The forks of the tines should be sitting over the wedge.

8. Lightly burnish the matrix against the adjacent tooth both buccally and lingually to further adapt it to form proper contours. **NOTE: Over burnishing may distort matrix and compromise contours.** The Pin tweezers has a built-in ball burnisher for this purpose. The matrix placement tab may be bent over the adjacent marginal ridge to facilitate restoration marginal ridge contour. Technique Tip: When restoring an MOD preparation, two rings may be placed facing opposite directions, one distally or both pointing in the same direction, usually mesially, one over the other placing the mesial matrix and ring first.

9. Placement of Restorative

Place restorative material as per restorative material manufacturer's instructions for use and establish proximal contact as per usual technique. Light cure and/or allow material to reach initial set. Gross contouring may be accomplished with matrix in place.

10. Matrix Removal

Following appropriate restorative curing/setting, remove ring(s) using forceps, remove wedge(s) and sectional matrices using pin tweezers. Finish contouring and verify appropriate proximal contact. Technique Tip: It is advisable to "post-cure" light cured restorative materials following matrix removal from all surfaces, including proximo-facial and lingual line angles.

Small parts under risk of ingestion/inhalation

Ensure matrix is grasped securely with pin tweezers, fingers or locking college pliers during withdrawal from the mouth. Recover all pieces from the mouth immediately.

4. HYGIENE AND DISPOSAL

Cross-contamination – To reduce risk of infection

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Reprocess reusable products as described below.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- During treatment clinics with patient contact should only handle items to be used.

NOTICE

Wrong cleaning or disinfection method.

Damage to components

- Ensure that phosphoric acid etchants are removed from the retaining ring tines prior to reprocessing
- Most ultrasonic cleaning solutions will not affect the Palodent® Plus rings, although chemicals in some cleaners may cause time discoloration. The discoloration is only cosmetic, but if that is a concern put the rings in a beaker with a water/detergent solution and place that in the ultrasonic cleaner to prevent contact with chemicals.
- Pre-Use Cleaning – Matrices, Wedges, Wedgeguards**
- Pre-use cleaning and disinfection may be omitted for single-use matrices, wedges and wedgeguards stored in their original packaging and protected from contamination.
- Proceed to packaging and steam autoclaving according to instructions below.

4.1 Palodent® Plus instruments and components

Process		✓ Acceptable
Cleaning	Mechanical (Ultrasonic)	✓
	Manual	✓
Disinfection	Automated (Washer/Disinfector)	✓
	Automated (Washer/Disinfector)	✓
Sterilization	Immersion	✓
	Steam Autoclave	✓
Single Use (Matrices, Wedges, Wedgeguards)	Dispose of after Use	✓
	Reusable/Not reprocessable	Dispose of if contaminated

Instructions for cleaning, disinfecting and sterilizing Palodent® Plus instruments and components

Forceps, Pin Tweezers, Retaining Rings, Matrices, Wedges, Wedgeguards		
Warnings	- High-level disinfection has not been validated as a terminal process for the Palodent® Plus Instruments and components. Disinfection is not recommended as a terminal process. Proceed to Sterilization following cleaning and any optional/intermediate disinfection process - Steam spore/occluding sterilization is appropriate and recommended for the Palodent® Plus Instruments and components. - Follow proper infection prevention activities, such as proper hand washing and donning new puncture and chemical gloves at pertinent steps. - Do not use hand wire brushes to clean the Palodent® Plus Instruments and components as scratching and mechanical degradation can occur. - Use only a disinfecting solution which is approved for its efficacy, EPA registered (and/or Health Canada approved) and use in accordance with the IFU of the disinfecting solution manufacturer - Do not allow the Palodent® Plus Instruments and components to exceed 137° C. - Do not use phenol-based glutaraldehyde - Always use a pH neutral instrument cleaning solution if authorized - Sterilize metal and plastic components in separate pouches to avoid damage to components - The integrity of the pouch has been compromised, the device should be processed again before use. - All devices should be cleaned and sterilized before each use. For the countries that require a three-step process, in which disinfection is required prior to sterilization, all devices must be cleaned, disinfected, and sterilized before use. - Repeated reprocessing has minimum effect on these instruments. End of life is normally determined by wear and damage due to use. - Cold liquid immersion disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use.	
Limitations on Reprocessing	- Do not allow residue or any form of contamination to dry on the device. Remove excess soil/gross debris by wiping, brushing, and/or rinsing under water. - Transport to processing area and follow validated processing instructions below. - It is recommended that the device be processed as soon as is reasonably practical following use.	
Initial Treatment at the Point of Use	- Do not attempt to disassemble. - Process the Palodent® Plus Instruments and components fully assembled. - The chemical manufacturer's instructions for the cleaning and disinfection solution should be followed, observing concentration rates, contact times, fill level, degassing processes, etc. - The recommended solution is Resurge® Instrument Cleaning Solution mixed to 0.5 oz/gal, at 35 - 40°C, FD 370 Cleaner, Dürr Dental AG, Germany or equivalent enzymatic cleaning solution - Immersion Palodent® Plus Instruments and components in an ultrasonic bath containing the prepared solution for at least 20 minutes unless otherwise specified by the time specified by the solution manufacturer - Rinse under tap water of at least drinking quality for a minimum of 30 seconds - Allow to air dry - Visually inspect for visible soil, and repeat the cleaning and disinfection procedure if necessary	

Manufactured for Dentsply Caulk 38 West Clarke Avenue Milford, DE 19963 USA 1-302-422-4511 www.dentsplysirona.com	Distributed by Dentsply Canada 525 Rowntree Dairy Rd Woodbridge, Ontario L4L 7K6 Canada
---	--

- Clean Palodent® Plus Instruments and components with a soft brush, tap water, and liquid soap until free of visible contamination. Pay special attention to device seams and interspersions.
- Remove cleaning solution residue with a cloth soaked with tap water or rinse under warm tap water for 2 minutes.
- Dry Palodent® Plus Instruments and components with a lint-free, single-use cloth or allow to air dry in a clean area.
- Visually inspect for visible soil, and repeat the cleaning procedure if necessary

- Disinfection by Wiping
- To disinfect Palodent® Plus Instruments and components, thoroughly wipe all surfaces, with a new single-use wipe in combination with an alcohol-based, tuberculocidal, quaternary ammonium cleaning and disinfecting solution
 - Use additional wipes as needed to ensure the entire outer surface, including seams and crevices, remain visibly wet for the contact time recommended by the solution manufacturer
 - Pay special attention to all seams and crevices while wiping
 - Rinse the Palodent® Plus Instruments and components under tepid running water (40 C, potable) for 30 seconds to remove disinfectant solution
 - Dry the Palodent® Plus Instruments and components using a dry, lint-free disposable cloth or disposable wipe

- A thermal disinfectant that complies with ISO 15883 may be used for automatic cleaning and disinfection.
- Use only properly maintained, calibrated, and approved washer disinfectors according to ISO 15883.
- Follow the manufacturer's instructions to select the solutions and programs
- Remove excess soil with a disposable cloth/paper wipe
- In case of gross contamination of Palodent® Plus Instruments and components, remove it using a soft brush and running water.
- Place Palodent® Plus Instruments and components in the washer-disinfector allowing water and detergent to enter into and drain out of device openings.
- Run washer-disinfector cycle with A0 value = 3000 (e.g. 5 min at ± 90°C) (e.g. Miele Vario TD) using appropriate detergents (e.g. Neodisher MediClean and Neodisher Z, both from Dr. Weigert, Germany or equivalent).
- After the cycle is complete, visually inspect the Palodent® Plus Instruments and components for degradation due to extended use and reprocessing. If the device is discolored, cracked, worn, distorted, etc., discard and do not use.
- Proceed to sterilization following Automated Cleaning and Disinfection.

- Drying
- The drying instructions are incorporated in the cleaning and disinfection section above.

- Maintenance, Inspection and Testing
- Visually inspect the Palodent® Plus Instruments and components for degradation due to extended use and reprocessing. If the device is discolored, cracked, worn, distorted, etc., discard and do not use.
 - Forceps, tweezers, and other devices require no additional maintenance and should not be oiled.
 - Ensure moving parts operate smoothly.
 - If a Palodent® Plus retaining ring is deformed by over stretching during placement or removal, the ring should be discarded.

- Packaging
- Ensure that the Palodent® Plus Instruments and components are completely dry before packaging for sterilization
 - A lint-free, disposable cloth may be used to dry Palodent® Plus Instruments and components
 - Place each Palodent® Plus device in a separate appropriately sized steam-sterilization pouch
 - Uses of FDA-cleared, ISO 11607 compliant Paper/ plastic steam sterilization pouches (e.g., AssurePlus® Sterilization Pouches) are recommended.

- Sterilization / Steam autoclaving
- The Palodent® Plus Instruments and components must be sterilized utilizing steam autoclave sterilization
 - Note:** Palodent® Plus Forceps, Pin Tweezers, Retaining Rings, Matrices, Wedges and Wedgeguards may be sterilized by Dynamic-Air-Removal/Pre-Vacuum steam autoclave.
 - The following Pre-Vacuum cycles may be used for Forceps and Pin Tweezers:
 - 134°C (273°F), for 3 minutes
 - 135°C (275°F) for 3 minutes

The following Pre-Vacuum cycles may be used for Retaining Rings, Matrices, Wedges and Wedgeguards:

- 132°C (270°F) for 4 minutes with a minimum drying time of 20 minutes
- 134°C (273°F) for 3 minutes with a minimum drying time of 18 minutes
- 135°C (275°F) for 3 minutes with a minimum drying time of 16 minutes

Note: Palodent® Plus Forceps, Pin Tweezers, Retaining Rings, Matrices, Wedges and Wedgeguards may be sterilized by Gravity Displacement steam autoclave.

- The following Gravity cycles may be used for Forceps and Pin Tweezers:
- 135°C for 10 minutes with a minimum drying time of 30 minutes.

- The following Gravity cycles may be used for Retaining Rings, Matrices, Wedges and Wedgeguards:
- 132°C for 15 minutes with a minimum drying time of 30 minutes.
 - 135°C for 10 minutes with a minimum drying time of 30 minutes.

An alternate method of sterilization is to place non-bagged instruments into the steam sterilizer and run at one of the above listed cycles

- Storage
- Do not remove from the pouch until ready for use to prevent contamination
 - Instruments sterilized non-bagged should be used immediately
 - Store at room temperature, away from moisture or excessive humidity.
 - To prevent contamination, store processed device in covered storage such as a drawer or cabinet until use.
 - Prior to use, inspect the pouch. If the integrity of the pouch has been compromised, the device should be processed again before use.
- Additional Information
- Prior to use, inspect the device. Discard any Palodent® Plus Instrument or component which has become damaged, worn, or distorted.
 - The instructions provided above have been validated by the manufacturer of the medical device as being capable of preparing a medical device for use/use. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and trained personnel in the processing facility, achieves the desired result. This requires verification and/or validation and routine monitoring of the process. Likewise, any deviation by the processor from the instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Use of any other reprocessing procedures/methods is at the discretion and sole responsibility of the practitioner.

- Manufacturer Contact
- Within the United States, call Dentsply Sirona at 1-302-422-4511. For areas outside the United States, contact your local Dentsply Sirona representative.

4.2. Disposal

- Dispose of in accordance with national and local regulations.
- During incineration, Wedges and WedgeGuards may produce irritating and highly toxic gases.
- Used matrices are contaminated and sharp. When handling or disposing, use precautions as when handling or disposing of other contaminated sharps.

5. LOT NUMBER AND CORRESPONDENCE

1. The following numbers should be quoted in all correspondence:

- Reorder Number
- Lot number on the package
- Copy of proof of purchase

2. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

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