ZOE B&T®

Base and Temporary Filling Material		
INSTRUCTIONS FOR USE	ENGLISH	2
Matériau de comblement temporaire et de base		
MODE D'EMPLOI	FRANÇAIS	4
Materiale di riempimento di base e temporaneo		
ISTRUZIONI PER L'USO	ITALIANO	6
Unterfüllungs- und temporäres Füllungsmaterial		
GEBRAUCHSANWEISUNG	DEUTSCH	8
Material de base y relleno temporal		
INSTRUCCIONES DE USO	ESPAÑOL	10

ZOE B&T®

Base and Temporary Filling Material

INSTRUCTIONS FOR USE - ENGLISH

Caution: This is a medical device. For dental use only by dental professionals.
USA: Rx only.

1. PRODUCT DESCRIPTION

ZOE B&T® Base and Temporary Filling Material is a polymer-reinforced zinc oxide-eugenol composition restorative material designed for intermediate restorations intended to remain in place for no longer than one year.

ZOE B&T® Base and Temporary Material may also be used as a base under restorative materials and cements that do not contain resin components.

When mixed properly, ZOE B&T® Material has a compressive strength in excess of 6300 psi (data on file).

ZOE B&T® Material complies with ISO 3107 Zinc Oxide/eugenol type II: for bases and temporary restorations.

1.1 Indications for use

- ZOE B&T® Material is indicated for:
- Intermediate, (less than 1 year) direct restorative filling material.
 - Base material under non-resin-containing permanent restorations.

1.2 Contraindications

- ZOE B&T® Material is contraindicated for use with patients who have a known hypersensitivity to eugenol.
- ZOE B&T® Material is contraindicated for use with patients who have a known hypersensitivity polymethylmethacrylate.
- ZOE B&T® Material is contraindicated for direct application to dental pulp tissue (direct pulp capping).
- ZOE B&T® Material is contraindicated as a base under resin-containing adhesives, restoratives or cements because eugenol may interfere with the hardening and/or cause softening of the polymeric (resin) components.

1.3 Delivery forms

- A powder/liquid version for hand mixing

1.4 Composition

Powder: Zinc oxide; Poly-Methyl Methacrylate (PMMA) powder; Zinc Acetate; Pigment
Liquid: Eugenol; Acetic acid

2. GENERAL SAFETY NOTES

Be aware of the following general safety notes and the special safety notes in other chapters of these Instructions for Use.

Safety alert symbol

This is the safety alert symbol. It is used to alert you to potential personal injury hazards. Obey all safety messages that follow this symbol to avoid possible injury.

2.1 Warnings

1. Liquid contains eugenol. Powder contains polymethacrylate. Components and mixed material may be irritating to skin, eyes and oral mucosa and may cause allergic contact dermatitis in susceptible individuals.
 - **Avoid Eye Contact** to prevent irritation and possible corneal damage. In case of eye contact, rinse immediately with plenty of water and seek medical attention.
 - **Avoid Skin Contact** to prevent irritation and possible allergic response. In case of contact, reddish rashes may be seen on the skin. If contact with skin occurs, immediately remove material with cotton and alcohol and wash thoroughly with soap and water. In case of skin sensitization or rash, discontinue use and seek medical attention.
 - **Avoid contact with oral soft tissues/ mucosa** to prevent inflammation. If accidental contact occurs, immediately remove material from the tissues. Flush mucosa with plenty of water after the restoration is completed and expectorate/evacuate the water. If inflammation of mucosa persists, seek medical attention.
2. Improved or reinforced zinc oxide-eugenol products such as ZOE B&T® Material are not suitable for use as root canal filling materials as inconsistent, incomplete or delayed resorption may occur.
3. ZOE B&T® Material may irritate pulp or soft tissue. Cover dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner) following manufacturer's instructions.
4. ZOE B&T® Material should not be used with patients who have a known hypersensitivity to any of the components.

2.2 Precautions

- This product is intended to be used as specifically outlined in the Instructions for Use. Any use of this product inconsistent with the Instructions for Use is at the discretion and sole responsibility of the practitioner.
- Insufficient data exist to support the use of ZOE B&T® Material in a root canal as a cement, dressing, medicament or repair material.
- Wear suitable protective eyewear, mask, clothing and gloves. Protective eyewear is recommended for patients.
- Contact with saliva, blood and sulcus fluid during application may cause failure of the restoration. Use adequate isolation such as rubber dam.
- Insufficient data exist to support the use of ZOE B&T® Material as a luting cement.
- Devices marked "single use" on the labeling are intended for single use only. Discard after use. Do not reuse in other patients in order to prevent cross-contamination.
- Always tightly close containers immediately after use.
- Increased temperature and/or humidity reduce available work time and accelerate set time. Material will set faster in the mouth than at the bench. Do not use bench cure rates to judge mouth cure rates.
- ZOE B&T® Material should be mixed in equal proportions (1 level scoop powder to 1 drop liquid, 6.0:1.0 by weight) for optimum performance. The mixed material should be homogeneous and streak free prior to application. Variations may affect the strength and durability of the material.
- Use in a well-ventilated area. Avoid inhaling dust particles from powder. (See Warnings)
- ZOE B&T® Material liquid bottle contains excess liquid for the amount of powder in one powder bottle. Providing excess allows for some natural liquid evaporation over repeated openings. Liquid bottle (containing remaining liquid) should be discarded and properly

disposed after all powder is used. Remaining liquid may become brown in color and should not be further used.

- Do not use instruments which have not been properly (re)processed. Before re-use, instruments must be properly cleaned/decontaminated to reduce the risk of cross-contamination.
- The bottles, scoop and dropper cannot be reprocessed. To prevent bottles, scoop and dropper from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that they are handled outside of the patient treatment operatory with clean/disinfected gloves. During treatment clinicians with patient contact should not handle bottles, scoop or dropper.
- Do not attempt to adjust material strength by incorporation of modifiers (e.g., petrolatum) or deviations from directed powder/liquid ratio.
- Product is not self-adhesive and may not be used in conjunction with conventional bonding techniques. Cavity preparation must provide for the mechanical retention of the product.
- Numerous factors can affect restoration longevity. Material is intended as an interim material. Check restoration marginal integrity and retention at every periodic examination.

Interactions:

- Do not use this product in conjunction with (meth)acrylate based dental materials since it may interfere with the hardening of the (meth)acrylate based dental materials.

2.3 Adverse reactions

1. Product may irritate the eyes and skin.
 - **Eye contact:** Irritation and possible corneal damage.
 - **Skin contact:** irritation or possible allergic response. Reddish rashes may be seen on the skin.
 - **Contact with Mucous Membranes:** inflammation. (See Warnings).
2. Product may cause pulpal effects. (See Contraindications).

2.4 Storage conditions

- Inadequate storage conditions may shorten the shelf life and may lead to malfunction of the product.
- Store with original (not dropper cap) caps tightly closed in a well ventilated place at temperatures between 10-24° C/50-75° F.
 - Allow material to reach room temperature prior to use.
 - Do not store with or near materials containing (meth)acrylate resin components.
 - Do not freeze.
 - Protect from moisture.
 - Do not use after expiration date.
 - Do not use if the lot number and/or expiration date are missing or illegible.

3. STEP-BY-STEP INSTRUCTIONS

3.1 Cavity preparation

- 3.1.1 Isolate the field of operation in preferred manner. Rubber dam is recommended.
- 3.1.2 Cavity preparation must provide for the mechanical retention of the product.
- 3.1.3 If caries is present, completely excavate using low-speed and/or hand instrumentation.
- 3.1.4 Place matrix and wedge where indicated.
- 3.1.5 Gently dry preparation with cotton pellet. Avoid desiccation.

3.2 Pulp protection

For direct and indirect pulp capping, cover the dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner, see complete Instructions for Use).

To reduce risk of pulp irritation

- Cover dentin close to the pulp (less than 1 mm) with a hard setting calcium hydroxide liner (Dycal® Liner) following manufacturer's instructions.

3.3 Dispensing ZOE B&T® Base and Temporary Filling Material

- 3.3.1 Before dispensing, fluff powder (tumble the powder bottle) to assure uniform bulk density of contents.
- 3.3.2 Using the dropper cap, dispense two drops (or more, for larger cavities) of liquid onto a clean, non-absorbent mixing pad or glass slab. After dispensing liquid, replace original cap immediately to prevent evaporation and subsequent contamination.
- 3.3.3 Dispense approximately 0.2 grams powder per drop dispensed onto the pad/slab.

NOTE: Keep powder and liquid separated until ready for spatulation.

NOTICE: Premature setting of material. Keep powder and liquid separated until ready for spatulation.

3.4 Mixing ZOE B&T® Base and Temporary Filling Material

- 3.4.1 Use a spatulation technique that quickly and thoroughly combines 50% of the dispensed powder into all the dispensed liquid.
- 3.4.2 Bring the remaining powder into the mix in 2 or 3 increments and spatulate thoroughly. The mix should be completed in approximately one (1) minute.
- 3.4.3 **Technique Tips:**
 - A proper mix will appear quite stiff and should be stropped (wiped vigorously) for 5 to 10 seconds.
 - A proper mix results in a working consistency that is smooth and adaptable. If a dry, stiff or rubbery mix is obtained, additional liquid may be incorporated. If a sticky and soft mix is obtained, more powder may be added. Slight increments of powder should be incorporated until proper consistency is reached.

3.5 Placement of mixed ZOE B&T® Base and Temporary Filling Material

Mixed material may be carried to the preparation with suitable placement instrument. Adapt, contour and shape with appropriate instruments.

Technique Tip: Dipping the instrument into dry powder left on the pad before using to adapt material will facilitate placement and minimize sticking.

3.5.1 Use as Intermediate Restorative Material

- Trim excess material. Material may be smoothed by lightly burnishing. Initial set occurs approximately 2 ½ – 3 ½ minutes from start of mix. Increased temperature, humidity, and powder/liquid ratio may accelerate set.
- Upon set, remove matrix band if present. Check and adjust occlusion as needed. Additional carving or occlusal adjustment may be accomplished with a round bur. Avoid gouging material with sharp burs at high speed. Additional polishing is not recommended.
- Periodically evaluate intermediate restoration until final restoration is indicated. The intermediate restoration should be removed and final restoration completed within 1 year of placement.
- 3.5.2 Use as Protective Base
- This product is a protective base for restorative material not containing polymeric resin components (e.g. final restoration of amalgam, gold)

Possible Interactions – To reduce risk of restoration failure

- The product should not be used in conjunction with resin-based adhesives, varnishes or restoratives since eugenol may interfere with hardening and cause softening of (meth)acrylate based dental materials.

- Apply mixed material to desired dentin surfaces and depth.
- Following completion of set (approximately 5 minutes from start of mix) cavity refinement may be completed.
- Complete restoration per restorative material manufacturer's instructions for use.

4. HYGIENE, PROCESSING AND DISPOSAL

Cross-contamination- To reduce risk of Infection

- Do not reuse single use products. Dispose of in accordance with local regulations.
- Bottles, scoop and droppers cannot be reprocessed. Dispose of contaminated bottles and/or droppers in accordance with local regulations.
- When not in use, all supplies and instruments should be in covered storage, such as drawers and cabinets, and away from potential contamination.
- To prevent bottles, scoop and droppers from exposure to spatter or spray of body fluids or contaminated hands it is mandatory that the bottles, scoop and droppers are handled with clean/disinfected gloves, outside of the patient treatment operator. Do not reuse bottles if contaminated.
- Dispensing materials with clean/disinfected gloves or hands or suitable transfer forceps in a separate room, bringing only what will be used into the operator is strongly recommended.
- During treatment clinicians with patient contact should only handle items to be used.
- To prevent bottles, scoop and droppers from exposure to spatter or spray of body fluids or contaminated hands, or oral tissues, use of a protective barrier is recommended. The use of protective barriers is an additional precautionary measure against gross debris but not against all contamination. The barrier is designed for single use and must be disposed of after each use in accordance with local regulations.

Incidental contact of the bottles, scoop and droppers with water, soap or a water-based hospital-level disinfection solution will not damage bottles, scoop and droppers. Do not allow any solution contact with contained material. Discard material that has been in contact with any fluid or non-sterile instrument.

Repeated liquid contact may damage label. Dry the bottles, scoop and droppers with a lint-free single-use cloth.

NOTE: Vigorous wiping can destroy the label. Wipe bottles gently.

4.2 Disposal

- Based on the information in the SDS, the waste generated must not be disposed of together with household garbage. Dispose of material as solid waste in a closed container.
- Dispose of contents/container to authorized hazardous or special waste collection point in accordance with any local regulation.
- Dispose of in accordance with national and local regulations.
- Contaminated containers and accessories that have been in contact with the patient must be disposed of as hazardous waste at risk of biological contamination.
- Dispose of disinfected mixing spatula, dispensing, and mixing vials and packages in normal waste stream in accordance with local regulations.
- Always wear gloves when handling packaging and accessories.
- See SDS for further details

5. LOT NUMBER, EXPIRATION DATE AND CORRESPONDENCE

1. Do not use after expiration date. Do not use if the lot number and/or expiration date are missing or illegible. ISO standard uses: "YYYY-MM-DD."
2. The following numbers should be quoted in all correspondence:
 - Reorder number
 - Lot number
 - Expiration date
3. Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.
4. Device Identification: Basic UDI-DI ++D002DENTALCEMENTSC6

©2025 Dentsply Sirona. All Rights Reserved. **507100WEB** (R 01/08/25)

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