Rinn® XCP-DS Fit®

Universal Sensor Holder

INSTRUCTIONS FOR USE	ENGLISH	2
Système de positionnement à un anneau et un bras		
MODE D'EMPLOI	FRANÇAIS	7
Sistema con un braccio e un anello di posizionamento		
ISTRUZIONI PER L'USO	ITALIANO	12
Ein-Ring- und Armpositionierungssystem		
GEBRAUCHSANWEISUNG	DEUTSCH	17
Sistema de posicionamiento con un anillo y brazo		
INSTRUCCIONES DE USO	ESPAÑOL	22

Instructions for Use - English

Indications for Use

XCP-DS FIT® holders are indicated for holding digital sensors for intraoral dental radiographs. XCP® and XCP-ORA® arms and rings are indicated for positioning and aligning intraoral dental radiographs.

Contraindications

None known.

Warnings

Sterilize device before use. Cold liquid disinfection/sterilization, chemical vapor sterilization, and dry heat sterilization methods have not been tested or validated for efficacy and are not recommended for use. Only sterilize with steam sterilization according to the provided sterilization instructions. DO NOT exceed 273°F (134°C). It is the user's responsibility to qualify alternative methods of disinfection and sterilization. Failure to follow properly validated sterilization processes and approved aseptic techniques may result in cross contamination. Any method of sterilization will shorten the life of plastic parts.

Precautions

Follow sterilization instructions carefully to prevent damage to plastic pieces during sterilization. Use sensor holders with a sensor cover. Do not pull on sensor cord during use.

Adverse Reactions

None Known.

Residual Risks

Failure to follow the validated reprocessing instructions provided may result in exposing the patient to cytotoxic residuals and/or cross-contamination.

Step-by-step Instructions



1. Place a protective cover over the sensor.
2. Load sensor into appropriate biteblock for the area to be imaged by holding sensor firmly and pressing edge of sensor against the movable clip.
3. Insert opposite sensor edge and release clip.
4. Insert positioning arm into the appropriate color-coded channel of the aiming ring and attach biteblock according to Positioning Guide.
5. Correctly assembled, the sensor should be centered when viewed through the aiming ring.
6. Position holder in patient's mouth and slide aiming ring close to patient's face.
7. Align x-ray tube to aiming ring.
8. Take exposure.
9. Change as needed for the next desired image.




Positioning Guide

Anterior




Always use a sensor cover to protect sensor.



Posterior

Bitewing

Endodontic

Infection Control Procedures

Reprocessing:

- Arm and Ring: Disassemble instrument(s) and separate all components prior to reprocessing. Plastic and metal components must be sterilized in separate pouches/packages.
- XCP-DS Fit Biteblocks: Separate the silicone band from the biteblock and retaining clip prongs, leaving the band loosely attached to the biteblock.
- The following procedures were validated for reprocessing always; observe all country specific requirements for reprocessing medical devices.

Cleaning: Manual

- Rinse the disassembled device with potable water while brushing with a soft-bristled brush to remove gross soil.
- Prepare an enzymatic instrument cleaning solution (such as Dentsply Sirona Resurge® Instrument Cleaning Liquid) and immerse the device per the manufacturer's instructions.
- Rinse the device with potable water for at least 30 seconds and dry using a disposable lint-free cloth.
- Visually inspect the device for visible soil. If soil is seen, repeat the cleaning procedure.

Cleaning: Automated - Mechanical (Ultrasonic Cleaners)

- Follow manufacturer instructions for instrument cleaning solution, observing concentration rates and contact times.
 - Note: Some instrument cleaning solutions also disinfect.
- May be cleaned using an ultrasonic water bath for up to 15 minutes.
- Prepare an enzymatic instrument cleaning solution per manufacturer's instructions paying attention to concentration rates and contact times. A pH neutral instrument cleaning solution such as Dentsply Sirona Resurge® Instrument Cleaning Liquid is recommended for ultrasonic bath.
- Adhere to the ultrasonic cleaning bath manufacturer's recommended liquid level.
- Place the device into the ultrasonic bath and cover for the indicated contact time; ensure the device is sufficiently covered.
- Remove the device from the ultrasonic cleaning bath and rinse with potable water for 30 seconds to remove cleaning agent.
- Inspect the device for visible debris. If debris remains, manually clean using a nylon brush wetted with the cleaning solution to ensure all visible contamination is removed. Rinse the device with potable water for 30 seconds to remove cleaning agent.
- Dry the device using a dry, disposable lint-free cloth or use a dental instrument dryer.

Cleaning and Disinfection: Automated – Instrument Washer-Disinfector

The washer-disinfector used must be approved by the manufacturer for the cleaning and disinfection of dental instruments and comply with ISO 15883-1/-2 (A0 value 600 or greater; e.g. 90 °C (194 °F) and 1 minutes holding time).

- Rinse the disassembled device with potable water while brushing with a soft-bristled brush to remove gross soil.
- Follow manufacturer instructions for instrument washer solution, observing concentration rates and contact times as well as the instrument washer's recommended cleaning and/or thermal disinfection cycle (1 minute at 90°C) with dry times..
- Visually inspect the devices for visible soil and repeat the cleaning and disinfection procedure if necessary.
- Washer disinfectors that conform to ISO 15883 have approved cleaning parameters in its program when the thermal disinfection cycle is selected. For reference, Rinn Holders were validated with the following programmed cleaning parameters:
 - 2 minutes pre-wash/rinse in cold tap water
 - 5 minutes enzymatic wash in warm water > 43°C
 - 2 minutes enzymatic neutralization in hot water at 60°C
 - 2 minutes rinse in hot tap
 - 10 minutes drying at 90°C
- Refer to the instructions for use of the washer-disinfector for further instructions.
- We recommend the use of an enzymatic cleaning agent, e.g. Valsure Enzymatic solution

Sterilization

Ensure that the device is thoroughly dry prior to sterilization. Enclose the device into an appropriately sized sterilization pouch/package intended for steam sterilization in accordance to ISO 11607. Place the pouched/packaged device into the sterilizer per the manufacturer's instructions. Sterilize the device according to one of the recommended methods listed below. Place bagged instruments into the steam sterilizer per manufacturer's instructions.

Sterilize metal and plastic components in separate pouches to avoid damage to the components. DO NOT exceed 273°F (134°C).

STEAM STERILIZATION CYCLES	STERILIZING TEMPERATURE	MINIMUM STERILIZING CYCLE TIME	MINIMUM DRY TIME
Pre-vacuum	≥ 270°F / 132°C	3 minutes	16 minutes
Steam-Flush Pressure-Pulse	≥ 270°F / 132°C	3 minutes	16 minutes
Gravity Steam Sterilization	≥ 270°F / 132°C	10 minutes	15 minutes

Warning: the device must be used immediately following sterilization if it is non-pouched/packed.

Storage


To maintain sterility, keep the device in the sterilization pouch/package until ready for use.

Any serious incident that has occurred in relation to this device should be reported to the manufacturer and the Competent Authority of the Member State.

Manufactured by:
Dentsply Rinn
1301 Smile Way
York, PA 17404 USA
www.dentsplysirona.com

See Website for Symbols Glossary: www.dentsplysirona.com/IFU

XCP-DS FIT® Kits

559900 XCP-DS FIT Hygiene Kit

559908 XCP-DS FIT Endodontic Kit

559909 XCP-DS FIT Hygiene Kit plus Endodontic Holders

559910 XCP-DS FIT Hygiene Kit with XCP® Arms and Rings

XCP-DS FIT® Anterior Holders

Fits Size 1, 1.5, 2

559901 Anterior Holder (2)

559921 Anterior Holder (10)

559906 Anterior - Thin Holder (2)

559926 Anterior - Thin Holder (10)

XCP-DS FIT® Posterior Holders

Fits Size 1, 1.5, 2

559902 Posterior Holder (2)

559922 Posterior Holder (10)

XCP-DS FIT® Bitewing Holders

Fits Size 1, 1.5, 2

VBW can also be used for HBW with Size 0

559904 Vertical Bitewing Holder (2)

559924 Vertical Bitewing Holder (10)

559903 Horizontal Bitewing Holder (2)

559923 Horizontal Bitewing Holder (10)

XCP-DS FIT® Endo Holders

Horizontal Fits Size 1.5, 2 / Vertical Fits Size 1, 1.5

559905 Endo Holder (2)

XCP-DS FIT® Size 0 Holders

Fits Size 0 H or V

559907 Size 0 Periapical Holder (2)

XCP-DS FIT® Replacement Bands

6 per pack

559913 Short (for Posterior, Vertical Bitewing, Endo and Size 0)

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