Conform Fit® Gutta-Percha Obturation Points

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FOR DENTAL USE ONLY

SINGLE USE

NON-STERILE – TO BE DISINFECTED BEFORE USE

INSTRUCTIONS FOR USE

REF. WGGPC, PTNGP, PTGGP, PTUGP, TSGP, VTBGP, PULG, TNGPF, RBGP, RMGP

1 PRODUCT DESCRIPTION

- Product description:
Conform Fit® Gutta-Percha Points are thermoplastic cones that are used to obturate the shaped, cleaned root canal. Each brand of Conform Fit® Gutta-Percha Points is designed to fit their corresponding file system.
- Intended purpose:
Gutta-Percha Points are intended to fill the shaped, cleaned and irrigated root canal space.
- Intended user:
Gutta-Percha Points are to be used only in a clinical or hospital environment, by qualified dental professionals such as general practitioners as well as endo specialists (endodontist) and dental assistants.
- Safety and clinical performance information and clinical benefits:
The clinical benefit for the patient being treated with Conform Fit® Gutta-Percha Points is to restore tooth integrity and function as part of endodontic treatment.

2 INDICATIONS FOR USE

- Application field:
Endodontic obturators.
- Indication for use:
Gutta-Percha Points are indicated for obturation of root canal.
- Intended patient population:
Patients who need root canal treatment.

3 CONTRAINDICATIONS

None known.

4 DELIVERY FORM

SKU	Product name	Packaging and content
WGGPC	WaveOne® Gold Conform Fit® Gutta-Percha	Box of 60 points
PTNGP	ProTaper Next® Conform Fit® Gutta-Percha	Box of 60 points
PTGGP	ProTaper Gold® Conform Fit® Gutta-Percha	Box of 60 points
PTGGPAST	ProTaper Gold® Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size F1 12 points of size F2 12 points of size F3 12 points of size F4 12 points of size F5
TSGP	TruShape® Conform Fit® Gutta-Percha	Box of 60 points
VTBGP	Vortex® and Vortex Blue® Conform Fit® Gutta-Percha	Box of 60 points
VTBGP04A2040	Vortex® and Vortex Blue® Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size 20/.04 12 points of size 25/.04 12 points of size 30/.04 12 points of size 35/.04 12 points of size 40/.04
VTBGP06A2040	Vortex® and Vortex Blue® Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size 20/.06 12 points of size 25/.06 12 points of size 30/.06 12 points of size 35/.06 12 points of size 40/.06
PULG	ProTaper Ultimate® Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size F1 24 points of size F2 24 points of size F3
TNGPF	TruNatomy® Conform Fit® Gutta-Percha	Box of 60 points
T00TNGPF00AST	TruNatomy® Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size S 36 points of size P 12 points of size M
T00TNGPF00AS4	TruNatomy® Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size S 24 points of size P 12 points of size M 12 points of size L
RBGP	RECIPROC® Blue Conform Fit® Gutta-Percha	Box of 60 points
TNSRBGPF00AST	RECIPROC® Blue Conform Fit® Gutta-Percha	Box of 60 points including: 36 points of size R25 12 points of size R40 12 points of size R50
RMGP	RECIPROC® Minima Conform Fit® Gutta-Percha	Box of 60 points

SKU	Product name	Packaging and content
TNSRMGPF00AST	RECIPROC® Minima Conform Fit® Gutta-Percha	Box of 60 points including: 12 points of size R20 24 points of size R25 12 points of size R35 12 points of size R45

5 COMPOSITION

- Components:
 - Zinc Oxide (45 - 55%)
 - Synthetic Gutta-Percha (20 - 30%)
 - Barium Sulfate (15 - 20%)
 - Titanium Dioxide (1 - 5%)
 - Yellow Colorant (< 2%)
 - Red Colorant (< 2%)
 - Remaining components (< 10%)

Physical and chemical properties

- Radiopacity: Equivalent to 6 mm of aluminum or more in thickness
- Melting Temperature: < 100°C (212°F)
- Nominal length: ≥ 28 mm (1.1 in)

6 COMPATIBLE DEVICES

Gutta-Percha Points are used in combination with:

- Spreaders and Pluggers.
- Shaping files and Paper Points.
- Calamus® and Gutta-Smart™.
- ThermaSeal Plus® Root Canal Sealer and Ribbon® Root Canal Sealer.
- AH Plus® Endodontic Sealer, ProRoot® Bio Sealer.

7 SAFETY NOTES

7.1 Warnings

- Gutta-Percha Points are single-use products. Single-use products shall neither be re-processed nor reused. Reusing these products increases the risk of cross-contamination and/or breakage.
- Use of the devices might lead to the following residual risks:
 - Dissatisfaction of the user and delay in treatment / patient discomfort.
 - Retreatment (initial infection never healed).

- Serious infection / serious health consequences for the patient (e.g. biological hazards such as infection, toxic or allergic reaction for the patient) have never been reported and are very unlikely but cannot be totally excluded.
- Fracture of the tooth due to obturation instruments.
- Patients with known history of allergy or hypersensitivity to the devices ingredients listed in this IFU should not use the device.

7.2 Precautions

- Safety and effectiveness of use have not been established in pregnant or breastfeeding women or in children.
- Inspect the packaging and product before use: if the packaging or the product is damaged, do not use the product and discard it.
- Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.
- Do not sterilize with heat.
- Do not use disinfecting solutions containing Phenol or any products which are not compatible with gutta-percha.
- Use a rubber dam system during the endodontic procedure.

7.3 Adverse reactions / Undesirable side effects

Temporary post-operative discomfort could arise generally as a result of receiving a root canal treatment.

7.4 Storage conditions / Handling

- Keep away from sunlight and heat.
- Maintain under room temperature ($\leq 30^{\circ}\text{C}$ (86°F)).
- Store in a dry area.

8 STEP BY STEP INSTRUCTIONS

The described Gutta-Percha Points are part of a complete system and should only be used in combination with the corresponding shaping files and paper points. They are identifiable with their system brand, size name and color code. The obturation step includes the use of Gutta-Percha Points and sealer. The Gutta-Percha Points are to be disinfected prior to use. After disinfection, other Gutta-Percha Points can be introduced into the canal depending on the technique used. A heated instrument can be used to soften the material and ensure perfect sealing and dimensional adaptation of the root canal. In the use of a heated instrument, follow the IFU of said instrument. Use an appropriate tool to cut the head, including color coding ring after the insertion of the Gutta-Percha Points in the root canal. The colored ring shall not remain in the root canal after insertion.

9 PROCESSING: HYGIENE, CLEANING, DISINFECTION AND STERILIZATION INSTRUCTIONS

- Immerse the Gutta-Percha Points in 5.25% NaOCl for 1 to 5 minutes at ambient temperature. Remove all the bubbles from the surface of the Gutta-Percha Points. Gently wipe the Gutta-Percha Points with a sterile gauze moistened with 70% isopropyl alcohol. Let the Gutta-Percha Points air dry.

- Do not use disinfecting solutions containing Phenol or any products which are not compatible with gutta-percha.
- Products are single use.

10 SAFE DISPOSAL OF THE DEVICE

Products shall be disposed of according to local regulations for the safe disposal of sharp and contaminated devices.

11 EXPIRATION DATE

Do not use after the expiration date indicated on the label.

12 SERIOUS INCIDENT REPORTING

Any serious incident in relation to the product should be reported to the manufacturer and the competent authority according to local regulations.

13 WARRANTY AND LIMITED REMEDY / LIMITATION OF LIABILITY

Disclaimer: The Instructions For Use have been validated by Dentsply Sirona. Users are solely responsible for any deviation from these instructions, and/or the use of alternative methods for endodontic treatment and/or processing when applicable. Dentsply Sirona accepts no liability for damage, injury, or any legal responsibility incurred directly or indirectly by the user due to a deviation from the Instructions For Use set forth above. The user shall observe safe and lawful practices including, but not limited to, those set forth in this document.

14 SUMMARY OF SAFETY AND CLINICAL PERFORMANCE

Information is available on IFU website dentsplysirona.com/ifu.

15 ADDITIONAL INFORMATION

See webpage dentsplysirona.com/ifu to:

- Get a free print copy of the IFU.
- Get an explanation of non-harmonized symbols for IFUs and labels.

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