NTI-Kahla GmbH Rotary Dental Instruments
Im Camisch 3 - D-07768 Kahla
Tel: +49 (0) 36424-573-0 Fax: +49 (0) 36424-573-29
nti@nti.de; www.nti.de

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 Reprocessing recommendation for the product families "Dental diamond grinding instruments" and "Diamond finishing strips" according to DIN EN ISO 17664-1

The following sections describe reprocessing before first use as well as after each use of the instruments of the product families "Dental diamond grinding instruments" and "Diamond finishing strips". This includes pre-treatment at the point of use, automated and manual reprocessing as well as steam sterilisation in a fractionated vacuum process.

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Procedures: The following reprocessing recommendation describes the procedures for automated and manual cleaning and disinfection as well as steam sterilisation in a fractionated vacuum process.

Products: Product families "Dental diamond grinding instruments" and "Diamond finishing strips"

1. Warnings	WARNING: All instruments must be cleaned, disinfected and sterilised before each use in accordance with the following reprocessing recommendations; this also applies in particular to first use after delivery. As part of your responsibility for the sterility of the instruments during use, ensure that only adequately validated methods specific for equipment and product are used for cleaning/disinfection and sterilisation. The validated parameters must be adhered to for each cycle. The devices and accessories used must be regularly maintained and checked. Proper and effective cleaning and disinfection is a prerequisite for effective sterilisation. Please observe the legal regulations applicable in your country as well as the hygiene regulations in the medical practice or hospital. This applies particularly to the different specifications for effective prion inactivation. Instruments must not be reused after opening the pulp and must be disposed of to avoid the transmission of prions. Isolate contaminated instruments from non-contaminated instruments until reprocessing, if necessary using suitable packaging and labelling. Personal protective equipment (PPE) must be worn when handling or working with contaminated or potentially contaminated materials, devices and equipment. Personal protective equipment must comply with the regulations of the respective country and must include at least the following: - protective clothing (impervious to fluids), - masks, - eye and face protection, - gloves and - overshoes. Exercise caution with instruments that have points or sharp cutting edges. Due to the intended purpose and the geometry/surface properties of the products, we recommend that they be sterilised under all circumstances in accordance with the classification "critical B" according to DIN EN ISO 17664-1. If unsuitable cleaning agents or disinfectants are used or if unsuitable procedures are employed, this can lead to the loss of warranty claims and have negative consequences for the instruments: - Damage - Corrosion - Discolouration of the product - the necessary hygiene not being ensured The ProxyStrip system must be disassembled into its individual parts (handle, ProxyStrip) for effective reprocessing.
--------------------	--

2. Limitations to reprocessing

2.1 Material resistance

Only use cleaning agents / disinfectants that are authorised for the reprocessing of dental instruments and observe the manufacturer's instructions. Do not use any fixing cleaning agents or disinfectants, otherwise blood contamination may become fixed.

When selecting the cleaning agent and disinfectant to be used, ensure that they do not contain the following ingredients:

- ▶ organic, mineral and oxidising acids
- ▶ strong alkalis (pH > 10.8) are not permitted, only neutral or weakly alkaline cleaning agents are recommended.
- ▶ alcohols, ethers, ketones and benzenes
- ▶ oxidants
- ▶ aldehydes

Do not clean the instruments with metal brushes or steel wool.

All instruments may only be exposed to temperatures up to a maximum of 138°C (280°F).

2.2 Reusability

After reprocessing the instruments with the maximum number of reprocessing cycles, no negative effects on the quality could be determined. The end of the service life is defined by the use of the instruments. Defective instruments, e.g. with missing diamond coating, fractured working parts, corroded surfaces and/or with detached colour markings as well as bent instruments, must be disposed of by the operator and must not be used.

The following number of reprocessing cycles has been validated for each product family:

Dental Diamond grinding instruments	up to 10 reprocessing cycles
Diamond finishing strips	up to 10 reprocessing cycles

Any further reprocessing and reuse as well as the use of damaged or contaminated instruments are the responsibility of the user. Any liability is excluded with failure to follow these instructions.

INSTRUCTIONS	
3. Pre-treatment at the point of use:	For pre-cleaning, use a plastic instrument bath (e.g. Fräsator) with an approved instrument disinfectant which, according to the manufacturer's instructions, does not contain any of the constituents stated in Section "2.1 Material resistance". Use separate containers for rotary instruments and diamond finishing strips. The disinfectant used must have tested efficacy (e.g. DGHM/VAH or FDA/EPA approval or CE marking), be suitable for instrument disinfection and be compatible with the instruments. For profiled instruments in particular, ensure that you have covered all surfaces and that there are no visible residues remaining and then rinse thoroughly under running water. If contamination is still visible, the instruments must be cleaned several times in accordance with the manufacturer's instructions to ensure complete removal. Instruments which cannot be cleaned of contamination must be sorted out at this point and disposed of professionally. NOTES: Pre-cleaning of the instruments should be performed immediately after treatment, taking personal protection into account. The aim is to prevent organic material and chemical residues from drying on the instruments and to avoid contamination of the environment. Please note that the disinfectant used for pre-treatment serves for personal protection and does not replace the subsequent disinfection step to be performed after completed cleaning. For manual removal of contamination do not use metal brushes or steel wool.
4. Storage and transport:	To prevent damage as well as contamination of the environment up to the point of further reprocessing, the instruments must be stored safely and transported in a closed container. Transport the instruments to the reprocessing room in a container provided for this purpose and place them in a basin with cleaning solution.
5. Preparation for cleaning:	All cleaning agents must be used at the dilution, temperature and contact time recommended by the manufacturer. Using the recommended temperatures, concentrations and contact time is important for optimum performance of the cleaning agents. NOTE: In accordance with the recommendation of the RKI, the cleaning bath must be freshly prepared at least every working day and changed immediately if visibly contaminated. In addition, the cleaning basin must be mechanically cleaned and disinfected every working day. Please also observe your country-specific requirements.
6. Automatic cleaning and disinfection	<u>6.1 Equipment:</u> When selecting the disinfectant, ensure that: - the disinfectant has proven efficacy (e.g. DGHM or FDA approval or CE marking in accordance with DIN EN ISO 15883) and - can run a tested programme for thermal disinfection (A0 value > 3000 or - for older appliances - at least 5 min at 90°C/ 194°F). The use of a washer-disinfectant compliant with EN ISO 15883-1 and -2 in combination with a suitable loading rack (e.g. corresponding small parts tray) is recommended. Follow the instructions for use issued by the manufacturer of the washer/disinfectant. NOTES: Load the washer/disinfectant (WD) according to the manufacturer's instructions.

	The concentrations, temperatures and exposure times specified by the manufacturer of the cleaning agent as well as the specifications for rinsing must be observed. Only use sterile or low-microbe (max. 10 microbes/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water (purified water/highly purified water). 6.2 Procedure: The instruments must be pre-cleaned (see Section "3 Pre-treatment at the point of use"). The instruments must be free of visible contamination prior to automated reprocessing. 1. Place the pre-cleaned instruments into the WD using a small parts tray. 2. Start the programme. 3. Remove the instruments from the WD immediately once the programme has ended. 4. Check whether post-drying is necessary and, if so, perform this in accordance with Section 8 "Drying". 5. Inspect and pack the instruments immediately after removal (see Section 9 "Maintenance, inspection and testing" and Section 10 "Packaging"). Reprocessing validation was performed under the following conditions in a Miele Professional PG8581 washer-disinfector using the neodisher® Mediclean Dental (0.5%, Dr. Weigert) detergent in the Vario TD programme: - Step 1 Pre-cleaning with cold tap water <40°C for 1 minute. - Step 2 Cleaning with 0.5% detergent 55°C for 5 minutes with demineralised water - Step 3 Rinsing with demineralised water for 1 minute - Step 4 Thermal disinfection with demineralised water >90°C for 5 minutes WARNING! Deviations from the validated disinfectant and the validated parameters or the described and validated procedure are the responsibility of the reprocessor/user. Never store the instruments in a wet condition after reprocessing. This can lead to corrosion and bacterial growth. If the instruments are not completely dry after completion of automated reprocessing, dry them manually (see Section 8 "Drying") and then store the instruments accordingly. NOTES: After thermal disinfection, a drying step can be included in the programme of the washer/disinfector. The validated parameters correspond to a process with an A0 value of >3000 s. It is recommended to only use procedures which meet or exceed this value.
7. Manual cleaning and disinfection:	7.1 Equipment: When selecting the cleaning agent and disinfectant to be used, ensure that the following preconditions are met: - Fundamental suitability for cleaning and disinfecting dental instruments - Suitability for ultrasonic cleaning (no foam formation)

- Use of disinfectants with tested efficacy (e.g. DGHM/VAH or FDA/EPA approval or CE marking)
- Compatibility of the chemicals used with the instruments (see Section 2 "Limitations to reprocessing")

NOTES:

Automatic cleaning / disinfection is preferred for reprocessing. Only use manual reprocessing if there is no option for automated reprocessing.

The prerequisite for using a combined cleaning agent/disinfectant is a very low initial load (no visible contamination) as the result of effective preliminary cleaning of the instruments.

It is essential that the concentrations, temperatures and exposure times stated by the manufacturer of the cleaning agent and disinfectant as well as the guidelines for final rinsing are complied with. Only use freshly prepared solutions, only sterile or low-microbe (max. 10 CFU/ml) as well as low-endotoxin (max. 0.25 endotoxin units/ml) water (purified water/highly purified water). Perform drying in accordance with Section 8 "Drying".

7.2 Procedure:

The instruments must be pre-cleaned (see Section "3 Pre-treatment at the point of use"). The instruments must be free of visible contamination **prior** to automated reprocessing.

Manual cleaning and disinfection were validated as being suitable for the reprocessing of instruments under the following conditions.

1. Prepare the cleaning/disinfection bath according to the instructions of the product manufacturer.
2. Place the pre-cleaned instruments into the cleaning bath (ultrasonic bath) so that the instruments are completely covered.
3. Then activate the ultrasound for the specified exposure time of the disinfectant used. The disinfectant neodisher Septo Plus (Dr. Weigert) with a concentration of 1.0% and an exposure time of 15 minutes in a 35 kHz ultrasonic bath type Bandelin Sonorex RK 514 H was used for validation.
4. Then remove the instruments from the cleaning bath and rinse them thoroughly at least five times (at least 1 minute) with low-microbe (max. 10 CFU/ml) and low-endotoxin (max. 0.25 endotoxin units/ml) water.
5. Dry the instruments according to Section 8 "Drying".
6. Check the instruments (see Section 9 "Maintenance, inspection and testing") for damage and visible residual contamination.
7. Pack the instruments as soon as possible after removal (see Section 10 "Packaging")

WARNING:

Deviations from the validated disinfectant and the validated parameters or the described and validated procedure are the responsibility of the reprocessor/user.

Never store the instruments in a wet condition after reprocessing. This can lead to corrosion and bacterial growth. If the instruments are not completely dry after completion of automated reprocessing, dry them manually (see Section 8 "Drying") and then store the instruments accordingly.

NOTE:

Chemical disinfection harbours the risk of disinfectant residues on the instruments and can affect the biocompatibility of the products, therefore thorough rinsing is absolutely essential.

8. Drying	Make sure to dry the instruments with medical compressed air (oil-free, low in microbes and particles) using a compressed air gun, air pressure: max. 10 bar, at room temperature or a clean, absorbent and lint-free cloth.						
9. Maintenance, inspection and testing	Check all instruments after cleaning resp. cleaning and disinfection (using a magnifying glass with 10-fold magnification) for corrosion, damaged surfaces, chipping, shape damage (e.g. instruments that are bent or no longer run smoothly), damage to the colour marking as well as contamination. Discard instruments that are damaged and instruments that cannot be successfully cleaned, disinfected and sterilised (for restrictions on reuse, see Section 2.2 "Reusability").						
10. Packaging	Pack the instruments in single-use sterilisation packaging (single or double packaging) which meets the following requirements: - Packaging suitable for steam sterilisation in accordance with EN ISO 11607 (temperature-resistant up to at least 138°C (280°F), sufficient steam permeability) - Adequate protection of the instruments against mechanical damage Select adequate size of the packaging to ensure the seal is not placed under tension. Also, do not use packaging that is too large to prevent the instruments from sliding around in the packaging.						
11. Sterilisation	The use of a steriliser in accordance with DIN EN 13060 Type B, EN ISO 17665 or EN 285 or ANSI AAMI ST79 is recommended. Sterilisation must be performed in packaging suitable for the sterilisation process. The packaging should comply with EN ISO 11607 (e.g. paper/laminate film). Moist heat / steam sterilisation in a minimum 3-fold fractionated vacuum process with sufficient product drying is the preferred and recommended sterilisation method. Sterilisation by steam sterilisation employing the gravitation method is not recommended. According to the following table, the drying time should not be less than 20 minutes. <table border="1"> <thead> <tr> <th>Country</th><th>Fractionated vacuum process</th></tr> </thead> <tbody> <tr> <td>EU</td><td>- 134°C (273°F) with at least 5 min holding time, 20 min drying time</td></tr> <tr> <td>Additional validated sterilisation conditions (observe the regulations of the respective country with regard to applicability)</td><td> - 132°C (270°F), min. 4 min holding time, 20 min drying time - 134°C (273°F), min. 3 min holding time, 20 min drying time - 135°C (275°F), min. 3 min holding time, 20 min drying time - 134°C (273°F), 18 min holding time, 20 min drying time* </td></tr> </tbody> </table> *The use of a prion programme is mandatory for sterilisation in several countries. Please observe the applicable national legal regulations. For fractionated pre-vacuum cycles, the country-specific requirements for the minimum parameter sets must be observed. WARNING: Do not use hot air sterilisation, radiation sterilisation, formaldehyde or ethylene oxide sterilisation, no chemiclaves as well as no plasma sterilisation. Never sterilise uncleaned instruments!	Country	Fractionated vacuum process	EU	- 134°C (273°F) with at least 5 min holding time, 20 min drying time	Additional validated sterilisation conditions (observe the regulations of the respective country with regard to applicability)	- 132°C (270°F), min. 4 min holding time, 20 min drying time - 134°C (273°F), min. 3 min holding time, 20 min drying time - 135°C (275°F), min. 3 min holding time, 20 min drying time - 134°C (273°F), 18 min holding time, 20 min drying time*
Country	Fractionated vacuum process						
EU	- 134°C (273°F) with at least 5 min holding time, 20 min drying time						
Additional validated sterilisation conditions (observe the regulations of the respective country with regard to applicability)	- 132°C (270°F), min. 4 min holding time, 20 min drying time - 134°C (273°F), min. 3 min holding time, 20 min drying time - 135°C (275°F), min. 3 min holding time, 20 min drying time - 134°C (273°F), 18 min holding time, 20 min drying time*						

	NOTE: The medical facility is responsible for the internal procedures for inspecting and packaging the instruments after they have been thoroughly cleaned to ensure vapour penetration and adequate drying. The user is responsible for validating the proper functioning of the steriliser.
12. Storage	Following sterilisation, the instruments must be stored in accordance with the instructions of the manufacturer of the sterile packaging.
13. Additional information	The above mentioned instructions have been validated by NTI-Kahla GmbH as being suitable for the preparation of "Dental diamond grinding instruments" as well as "Diamond finishing strips" for reuse. It is the responsibility of the medical facility to ensure that reprocessing is performed with the appropriate equipment and materials and that the staff in the reprocessing unit has been adequately trained to achieve the desired result. This necessitates validation and routine monitoring of the process. Furthermore, any deviation by the reprocessor from the recommendations should be properly evaluated for effectiveness and potential adverse consequences. Users must establish an appropriate cleaning protocol for the reusable medical devices used at their sites in accordance with the recommendations of the device manufacturer and the detergent/disinfectant manufacturer. Due to the numerous variables involved in sterilisation/decontamination, each medical facility should validate the sterilisation/decontamination process (e.g. temperatures, times) used with their devices.
14. Contact with the manufacturer	NTI-Kahla GmbH Rotary Dental Instruments Im Camisch 3 07768 Kahla Germany Tel.: +49(0)36 424-573-0 www.nti.de

 NTI-Kahla GmbH Rotary Dental Instruments
Im Camisch 3 - D-07768 Kahla
Tél. : +49 (0) 36424-573-0 Fax : +49 (0) 36424-573-29
nti@nti.de; www.nti.de

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