

TECHNICAL DATA SHEET
5 % Sodium Fluoride White Varnish

Intended Use

For the prevention of caries.

Active Ingredients

Sodium Fluoride 5.0% w/w (2.26% fluoride ion equivalent)

Characteristic	Specification
Description	Plastic tray with peel-off foil lid containing one brush and a gel preparation
Appearance	Viscous liquid, uniform consistency
Color	Off-white, slightly opaque
Format / Size	0.4 ml/Tray
Treatment	Non-Sterile

Additional Information

Storage Condition	Shelf Life	Country of Origin
Store under normal storage condition (dry, well ventilated, premises at 15 - 28 °C). Keep away from open flame and direct sunlight.	2 years	USA

Appendix 1 List of Item Codes and Attributes

Appendix 1 (List of items)

AMD Medicom Inc

Item Code	Description	Packaging Information	Specific Attribute
1031-C200	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	200 Trays/Box Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Caramel
1031-C30	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	30 Trays/Box, 6 Boxes/Pack, 180 Trays/Pack Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Caramel
1031-M200	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	200 Trays/Box Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Mint
1031-M30	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	30 Trays/Box, 6 Boxes/Pack, 180 Trays/Pack Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Mint
1031-ST200	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	200 Trays/Box Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Strawberry
1031-ST30	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	30 Trays/Box, 6 Boxes/Pack, 180 Trays/Pack Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Strawberry
1031-CH30	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	30 Trays/Box, 6 Boxes/Pack, 180 Trays/Pack Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Cherry
1031-CH200	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	200 Trays/Box Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Cherry
1031-BG30	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	30 Trays/Box, 6 Boxes/Pack, 180 Trays/Pack Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Bubblegum
1031-BG200	Medicom Duraflor Ultra 5% Sodium Fluoride White Varnish	200 Trays/Box Outer Packaging: plastic shrink-wrap Each box contains 1 Pad and 1 Product Insert and Instructions for use	Flavor: Bubblegum