SciCan

OPTIM 33TB

One-Step Cleaner & Intermediate Disinfectant

Active Ingredient: **0.5% Hydrogen Peroxide**



1 MINUTE
TUBERCULOCIDAL
CONTACT TIME

The Perfect Balance

No compromise between safety and efficacy.

OPTIM 33TB is one of the fastest and most effective cleaner & intermediate level disinfectants. OPTIM 33TB is safe; it does not require any GHS hazard pictograms, signal words or statements. OPTIM 33TB provides you with confidence that you are compliant and protecting your patients, staff, equipment and the environment; the perfect balance.

Effective Cleaner

- Removes organic debris while disinfecting, saving you time and money[1]

Fast Broad-Spectrum Disinfectant

- 1 minute tuberculocidal, bactericidal, virucidal contact time
- 3 minute fungicidal contact time
- 1-minute contact time against COVID-19
- Demonstrated 1-minute virucidal efficacy extends to new and emerging pathogens
- Ensures compliance, saving you time and money[2]

Safe for yourself and your patients

- Non-irritating to skin, eyes and respiratory tract
- Non-toxic to reproduction; child development and fertility
- No VOCs (volatile organic compounds)

Environmentally-Friendly

- Active ingredient, Hydrogen Peroxide, breaks down into water and oxygen
- No active residues, reducing environmental impact

Time and Cost Efficient

- Half as many wipes as the leading competitor are needed to keep your surface wet long enough to achieve disinfection[3]

Compatible

- Designed to be compatible with a wide range of materials[4]

Made in Canada
Fabriqué au Canada

COLTENE

Fast Contact Times

A contact time is the amount of time a surface must stay wet for disinfection to occur. OPTIM 33TB's contact times are faster which ensures compliance and provides confidence that disinfection has been achieved.

OPTIM 33TB has been proven effective against difficult to kill pathogens such as Norovirus, giving you confidence that your facility is protected against pathogens of concern.

Surface Disinfectant Contact Time*

Product Name	Active Ingredient	Mycobacteria	Non-enveloped Viruses				Enveloped Viruses					Bacteria				Fungi
		TB	Adeno	Polio	Noro	Rota	RSV	HIV	HBV (Hep.B)	HCV (Hep.C)	Influenza A	MRSA	VRE	E.coli	PSA	
OPT*IM* 33TB	Hydrogen Peroxide	1		1	1	1		1	1	1	1	1	1	1	1	3
Birex SE	Phenolics	10	10					10			10	10		10	10	10
Patterson pdCare	Quats-alcohol	3						2	2	2	2	2	2		3	3
Micro-Kleen 3	Quats-alcohol	3		3	3	3		3	3	3	3	3	3		3	3
CaviWipes	Quats-alcohol	3						3		3	3	3	3	3	3	3
CaviWipes 2.0	Quats-alcohol	2	2		2	2		2	2	2	2	2	2	2	2	2
CaviWipes HP	Hydrogen Peroxide	1	1		1	1		1	1	1	1	1	1	1	1	1
Sani-Cloth Super	Quats-alcohol	2	2			2	2	2	2	2	2	2	2	2	2	2
Germicide 3	Quats-alcohol	3		3				3	3	3	3	3	3		3	3
BioText	Alcohol Anhydrous	5													3	3
BioSurf	Alcohol Anhydrous	1													20 sec	10 sec
tb Minuteman	Alcohol Anhydrous	1		1				1			1	1	1		1	1
Clorox HP Wipes	Hydrogen Peroxide	1			3	1	30 sec	30 sec	30 sec	30 sec	1	1	30 sec	30 sec	30 sec	5

* Contact time in minutes. If blank, no published information was available. Information was obtained from publicly available material in April of 2024. Information will be updated if there are changes in claims.



Ordering Information for OPT*IM* 33TB

Ready-to-use Pre-saturated Wipes	
SCI33W12	OPTIM 33TB RTU Disinfectant Wipes, Unscented (12 x 160 sheets)
SCI33W12-ALOE	OPTIM 33TB RTU Disinfectant Wipes, Aloe Scent (12 x 160 sheets)
Liquid	
SCI334X4	OPTIM 33TB RTU Disinfectant, Unscented, 4L Bottle (4x4L)
SCI334X4-ALOE	OPTIM 33TB RTU Disinfectant, Aloe Scent, 4L Bottle (4x4L)
Dry Wipes Refill & Accessories	
SCI33R18	Dry Wipes Refill, 6" x 7" (18 x 160 sheets / case)
SCI33R18XL	Dry Wipes Refill, Large, 10" x 10" (18 x 60 sheets / case)
FSA01	Sprayer for 1L Bottles

For more information contact:
SciCan Ltd.
1440 Don Mills Road, Toronto, Ontario M3B 3P9 Canada
Tel: 416-445-1600 Fax: 416.445.2727
Toll Free: 1.800.667.7733
customerservice@scican.com

[1] John A. Molinari, Ph.D., "OPTIM 33TB Cleaning Study," Ann Arbor: Dental Advisor, 2010.
[2] John A Molinari, "Environmental Surface Wetness Test: Comparison of Disinfectant Wipes," Ann Arbor: Dental Advisor, 2015.
[3] John A Molinari, "Environmental Surface Wetness Test: Comparison of Disinfectant Wipes," Ann Arbor: Dental Advisor, 2015.
[4] John Molinari, Ph.D., Peri Nelson, B.S. Effect of OPTIM 33TB on Tensile Properties of Dental Chair Fabrics. Ann Arbor: Dental Advisor, 2016.

www.scican.com

OPTIM is a registered trademark of SciCan Ltd.

Manufactured by: Virox Technologies Inc. | 2770 Coventry Rd, Oakville, ON L6H 6R1 | Phone: (905) 813-0110

COLTENE

SD-248 CA EN R9